# Latex Gloves Powdered

**Sizes** XS - XL



## **Glove Specifications**

Size	Weight g/pc	Length mm	Palm Width mm	Cuff	Thickness Palm	Finger
XS	3.8 ± 0.2	240	75 ± 10	0.07 ± 0.02	0.085 ± 0.02	0.10 ± 0.02
S	4.3 ± 0.2	240	80 ± 10	0.07 ± 0.02	0.085 ± 0.02	0.10 ± 0.02
М	4.8 ± 0.2	240	95 ± 10	0.07 ± 0.02	0.085 ± 0.02	0.10 ± 0.02
L	5.3 ± 0.2	240	111 ± 10	0.07 ± 0.02	0.085 ± 0.02	0.10 ± 0.02
XL	5.8 ± 0.2	240	115 ± 10	0.07 ± 0.02	0.085 ± 0.02	0.10 ± 0.02
Tensile strength before ageing (all gloves)				Load @ peak	MPA	Min.14
Elongation before ageing (all gloves)				600%		



AQL4

EN445 Parts 1, 2 & 3

**EN384** 

Class 1 Medical Device

MDD 93/42/EEC

#### **Performance Requirements**

Test	Related Defects	Inspection Level	AQL
Water Leak at Final Release	Holes	G-1 Single Sampling Plan	4

### **Manufacturing Standards**

Manufacturing Standards	Incompliance with ISO 9001:2008
Presentation	Standard: 100 gloves per dispenser box;
	10 dispenser boxes per Carton
Glove Care	Should not store under direct sunlight or at high heat and humid conditions. Single use only.



#### **General Details**

Material	Natural Rubber Latex		
Surface Finish	Powdered		
Colour	Natural		
Geometry	Ambidextrous		
Texture	Fully Textured (Palm Textured)		
Cuff and Finishing	Beaded		